

# FDA Submission

**Your Name:** Ahmed M. Fathy

**Name of your Device:** Pneumonia Detection from Chest X-Rays software.

## Algorithm Description

### 1. General Information

#### Intended Use Statement:

- This algorithm is intended for use on both women and men from the ages of 5-85 who have been administered a chest X-Rays pneumonia study using PA or AP position on a X-Rays machine.

#### Indications for Use:

- for assisting radiologists in the **confirming** of pneumonia on chest X-Rays images.


#### Device Limitations:

- The Algorithm has a high false negative rate, so it is not recommended for using in screening or prioritization.

#### Clinical Impact of Performance:

- This algorithm has a high precision (low false positive rate), so it can be used for helping radiologists in **confirming** pneumonia diagnosis.

### 2. Algorithm Design and Function



**DICOM Checking Steps:**

- Patient's age between 5 and 85.
- Patient's examined body part is chest.
- Patient's position is PA or AP

**Preprocessing Steps:**

- Images are scaled by 1/255
- Images are resized to [224,244,3]

**CNN Architecture:**

Layer type	Output Shape	Param	note
InputLayer	(None, 224, 224, 3)	0	VGG16(freezed)
Conv2D	(None, 224, 224, 64)	1792	VGG16(freezed)
Conv2D	(None, 224, 224, 64)	36928	VGG16(freezed)
MaxPooling2D	(None, 112, 112, 64)	0	VGG16(freezed)
Conv2D	(None, 112, 112, 128)	73856	VGG16(freezed)
Conv2D	(None, 112, 112, 128)	147584	VGG16(freezed)
MaxPooling2D	(None, 56, 56, 128)	0	VGG16(freezed)
Conv2D	(None, 56, 56, 256)	295168	VGG16(freezed)
Conv2D	(None, 56, 56, 256)	590080	VGG16(freezed)
Conv2D	(None, 56, 56, 256)	590080	VGG16(freezed)
MaxPooling2D	(None, 28, 28, 256)	0	VGG16(freezed)
Conv2D	(None, 28, 28, 512)	1180160	VGG16(freezed)
Conv2D	(None, 28, 28, 512)	2359808	VGG16(freezed)
Conv2D	(None, 28, 28, 512)	2359808	VGG16(freezed)
MaxPooling2D	(None, 14, 14, 512)	0	VGG16(freezed)
Conv2D	(None, 14, 14, 512)	2359808	VGG16(Tunned)
Conv2D	(None, 14, 14, 512)	2359808	VGG16(Tunned)
Conv2D	(None, 14, 14, 512)	2359808	VGG16(Tunned)
MaxPooling2D	(None, 7, 7, 512)	0	added

Flatten	(None, 25088)	0	added
Dense	(None, 1024)	25691136	added
Dropout	(None, 1024)	0	added
Dense	(None, 512)	524800	added
Dropout	(None, 512)	0	added
Dense	(None, 2)	1026	added

### 3. Algorithm Training

#### Parameters:

- Adam optimizer was used with learning rate 0.0001 and decay 1e-5.
- Binary cross entropy was used for the loss.
- Accuracy matrix was used.
- Number of epochs 15.

#### Types of augmentation used during training:

- rescale=1. / 255.0
- horizontal\_flip
- Height\_shift\_range = 0.1
- Width\_shift\_range = 0.1
- Rotation\_range = 20.0
- Shear\_range = 0.1
- Zoom\_range = 0.1

#### Batch size:

- Training = 16
- Validation = 32

#### Optimizer learning rate:

- 0.0001

#### Layers of pre-existing architecture that were frozen:



- The first 15 layers In VGG16 model (see the CNN Architecture).

#### Layers of pre-existing architecture that were fine-tuned:

- The 3 convolution layers after 15 layer and the max pooling layer (see the CNN Architecture).

#### Layers added to pre-existing architecture:



- Flatten - Dense - Dropout - Dense - Dropout - Dense (see the CNN Architecture).






### Threshold and Explanation:

The threshold was chosen to be **0.51**. This achieves **0.8 precision** and **0.44 recall** and makes the algorithm efficient in **confirming** the pneumonia cases.

#### 4. Databases







### Description of Training Dataset:

- The training dataset is 80% of the total dataset. It was balanced and has 2290 cases.

### Description of Validation Dataset:

- The validation dataset is 20% of the total dataset. It was balanced and has 572 cases.

## 5. Ground Truth

- The disease labels were created using Natural Language Processing (NLP) to mine the associated radiological reports.
- The biggest limitation of this dataset is that image labels were NLP-extracted so there could be some erroneous labels but the NLP labeling accuracy is estimated to be >90%.

## 6. FDA Validation Plan

### Patient Population Description for FDA Validation Dataset:

- To validate the algorithms, I would collect a validation set that was made up of 2d X-rays chest images with both PA and AP positions for both women and men between the ages of 5 and 85. I would also want to make sure that the distribution of pneumonia in my validation set was reflective of the distribution of the density that is seen in the real world.

### Ground Truth Acquisition Methodology:

- As the diagnosis of pneumonia from chest X-rays is difficult, **the silver standard approach** of using several radiologists was used.

### Algorithm Performance Standard:

- The algorithm F1 score is 0.57, which is equal to the mean of F1 scores published in [this work](#).